

Positions you well.

GMFCS
Level IV-V

kg
max. 100



75-185 cm

0°-90°



I am your height adjustable supine standing frame.

With me, your transfer from the lying to the standing position is safe and comfortable – even if you are not actively involved. With my various adjustment options, I adapt very precisely to your specific body shape – from head to toe. Regular standing has a positive effect on your body structures and functions and promotes your activities. In doing so, I support you and at the same time offer you security. With me you have a new perspective – you can be at eye-level with your friends and can be moved around. This makes participation in different situations and environments during standing possible and enjoyable.

Height adjustment

is easily done via motor or with gas pressure spring support.



Growth adjustments

can easily be made by adjusting the length of the back and legs.







Push bar
for easy transfer.

Headrest
three-dimensionally
adjustable.

Upper body support
offers comfortable support.

Trunk and pelvic support
supports you and facilitates
a stable and upright
standing position.

Therapy tray
with padding offers plenty
of space for activities

Knee angle
adjustable up to 65°.

Flexible knee positioners (FKP)
are individually three-dimensionally
adjustable. Length-adjustable
contact straps adapt optimally to
the shape of the knee.

Footrest
individually height adjustable with adjustment
for dorsi/plantar flexion and eversion/inversion,
with heel edges incl. strap with magnetic closure.

5 Castors
secure, stable and easy to move



You can also easily adjust my height and angle. So nothing stands in the way of playing with my friends.



Precise tilting of the lying surface is possible from 0-90° via gas pressure spring.





Headrest with adjustable side guide

is flexible, adaptable, practical and offers great therapeutic benefits.



Flexible knee positioners (FKP)

are individually three-dimensionally adjustable. Length-adjustable contact straps adapt optimally to the shape of the knee.

An optimal foot position is achieved by individually adjustable 3D footrests. These can be lowered to the floor.



Perfect leg guidance is made possible by legs adjustable in length and angle.





Precise abduction of the legs
is possible from 0-30°
via adjusting screw.

Data		Size 1	Size 2	Size 3	Size 4
Body size		75 - 115 cm	100 - 135 cm	120 - 160 cm	145 - 185 cm
Back	Hip/thorax width	14 - 27 cm	17 - 31 cm	20 - 37 cm	27 - 46 cm
	Length of lying surface	30,5 - 43 cm	38 - 48 cm	49 - 61 cm	58 - 69 cm
Legs	Sole to knee	20,5 - 30,5 cm	27,5 - 37 cm	32,5 - 45 cm	42 - 54 cm
	Knee to crotch	12 - 18,5 cm	16,5 - 24 cm	20 - 29,5 cm	24,5 - 34 cm
	Knee angle	0° - 65°			
	Hip angle	0° - 45°			
	Abduction	0° - 30°			
Level compensation		0 - 5 cm			
Feet	Equinus	0° - 20°			
	Flat foot	0° - 15°			
	Rotation	0° - 15° per direction			
Lying surface height (for transfer)		53,5 - 73,5 cm	60 - 85 cm	71,5 - 93,5 cm	82,5 - 104,5 cm
Tilt		0° - 90°			
Transport dimensions (L x W x H)		95 x 51 x 75 cm	105 x 61 x 85 cm	125 x 76 x 97 cm	152 x 77 x 110 cm
max. load		30 kg	60 kg	80 kg	100 kg
Weight		30 kg	35 kg	80 kg	90 kg
HMV-No.		28.2902.1057			



You can find more information on accessories
in the price list.



Scan this QR code
and watch videos about
the product directly on our
homepage.

