

by Rifton



## **TECHNICAL SPECIFICATIONS**

	Small Wave	Medium Wave	Large Wave
User's Overall Height	(76 – 117 cm)	(102 – 142 cm)	(127 – 188)
Seat Height	2 – 7	2 – 7	2 – 7
	With tub stand: 12.5 – 18.5	With tub stand: 12.5 – 18.5	With tub stand: 12.5 – 18.5
	With shower stand: 27.5	With shower stand: 27.5	With shower stand: 27.5
	With tub transfer base (with standard leg inserts): 21 – 26	With tub transfer base (with standard leg inserts): 21 – 26	With tub transfer base (with standard leg inserts): 21 – 26
	With tub transfer base (with long leg inserts): 24.5 – 29.5	With tub transfer base (with long leg inserts): 24.5 – 29.5	With tub transfer base (with long leg inserts): 24.5 – 29.5
	With tub transfer base (with long leg inserts): 24.5 – 29.5	With tub transfer base (with long leg inserts): 24.5 – 29.5	With tub transfer base (with long leg inserts): 24.5 – 29.5
Seat Depth	10.5	13	17
Seat Width	15	15	17.5
Back Height	23	26.5	32
Length of Mesh Area when flat	35 (without calf rest)	41 (without calf rest)	50.5 (without calf rest)
Length of Mesh Area when Flat	42½ (with calf rest)	49 (with calf rest)	61 (with calf rest)
Base Width	18	18	21
Base Length	legs down 24 / legs flat 29.5	legs down 24 / legs flat 29.5	legs down 24 / legs flat 29.5
Maximum Working Load (lbs)	75	150	250
Back and Seat Angle Adjustments	0°, 15°, 30°, 45°, 60°, 75°	0°, 15°, 30°, 45°, 60°, 75°	0°, 15°, 30°, 45°, 60°, 75°



Glencar Enquire about this product



## The Wave

by Rifton



## **TECHNICAL SPECIFICATIONS**

Tub Stand Dimensions		
Height (legs flat)	8	
Length (legs flat)	39	
Length (legs down)	29.5	
Width	17.5	

Shower Stand Dimensions		
Height	22	
Length	34	
Width	27	

Tub Transfer Base Dimensions			
Max. Tub Height for Transfer Base with Standard Leg Inserts	20		
Max. Tub Height for Transfer Base with Long Leg Inserts	24		
Length	30		
Width	28		



Glencar Enquire about this product