

The HTS Toileting and Showering System

by Rifton



TECHNICAL SPECIFICATIONS

		Size 1 /Small	Size 1 / Medium	Size 3 / Large
Width	Mobile Base	53 cm	53 cm	57 cm
	Stationary Base	51 cm	51 cm	53 cm
	Hi-Lo Base	67 cm	67 cm	71 cm
Seat Height	Mobile Base	39 – 55 cm	39 – 55 cm	44 – 58 cm
	Stationary Base	33 – 43 cm	36 – 44 cm	38 – 48 cm
	Hi-Lo Base	51 – 78 cm	51 – 78 cm	51 – 78 cm
	Seat width	20* / 25* / 30 cm	25* / 30* / 36 cm	33* / 38* / 43 cm
	Seat depth	20 / 23 / 25 / 28 cm	29 / 32 / 34 / 37 cm	37 / 40 / 42 / 44 cm
	Armrest height	13,5 cm	17 cm	20 cm
	Back width	29 cm	34 cm	41 cm
	Lower leg length, toilet seat	22 – 34 cm	27 – 39 cm	34 – 44 cm
	Back height	42 cm	50 cm	62 cm

The HTS Toileting and Showering System

by Rifton



TECHNICAL SPECIFICATIONS

	Size 0 Nursery	Size 0	Size 1 Nursery	Size 1
Max. Toilet Height over which the Product can be Positioned	Mobile Base	52 cm	52 cm	52.5 cm
	Stationary Base	38 cm	38 cm	42 cm
Length	Mobile Base	74 cm	74 cm	81 cm
	Stationary Base	71 cm	71 cm	74 cm
	Hi-Lo Base	84 cm	84 cm	87 cm
	Front edge of hygiene opening to rear edge of mobile Base	48 cm	48 cm	48 cm
	Back angle	(-)6°, 6°, 18°	(-)6°, 6°, 18°	(-)6°, 6°, 18°
	Seat angle	(-)5° – 25°	(-)5° – 25°	(-)5° – 25°
		(-)15° – 15°	(-)15° – 15°	(-)15° – 15°
	max. load	34 kg	68 kg	113 kg

* can be adjusted using optional thigh guide